

RESPONDER[®] Latex-Free CPR Masks

AEDSuperstore[®]

An ALLIED 100[™] Company | World's Largest Automated External Defibrillator Source



Reliable, Affordable Protection for Rescuers

Our compact, precision-designed mask meets or exceeds the quality of the name brands at a low value price.

Order Today! Ships Today!

Form fitting heart shape adapts to both adult and child proportions. Pointed side up on adults. Pointed side down on children.

Clear, see-through face piece affords full view of the victim.

Disposable, one-way valve protects the responder by reducing the risk of contact with contaminants. Can also be used with training valves.

Oxygen barb for the administration of supplemental oxygen, if desired.

Elastic strap helps keep the mask in place. Frees you to do chest compressions.

BVM Compatible. With valve removed, fits onto bag valve mask resuscitators.

Soft, pliable cuff forms excellent seal in all weather conditions.



Ships collapsed.



Pop up to use.



AMP 0826
Responder CPR Mask
in Hard Case



AMP 0828
Responder CPR Mask
with Gloves and Wipe
in Hard Case



AMP 0823
Responder CPR Mask
with Gloves and Wipe
in Soft Case



AMP 0825
Responder CPR Mask
in Poly Bag

AEDSuperstore.com

Phone 800-544-0048

Email help@aed.com