

# CARDIAC SCIENCE AEDs

**G3** *third generation*

## FIRSTSAVE<sup>®</sup> AED G3 automated external defibrillator



### INTRODUCING THE CARDIAC SCIENCE FIRSTSAVE AED G3 TECHNOLOGY AND PERFORMANCE ... IN A VALUE PRICED DESIGN

The next generation FirstSave AED has arrived. Technology and value – the FirstSave AED G3 offers customers our leading edge AED technology in a value priced AED design. The new Cardiac Science FirstSave AED G3 features:

- Simple, easy to use, one-button operation
- Patented technology performs daily, weekly, and monthly self-tests
- More instructive voice prompts guide user through rescue
- Outstanding voice prompt quality and clarity
- Lightweight in a new user-friendly form factor

The FirstSave AED G3 also incorporates our advanced technology features such as the Company's patented RHYTHMx<sup>®</sup> analysis software and STAR<sup>®</sup> biphasic defibrillation energy waveform. Industry leading technology and performance features in a value priced, user-friendly form factor leads the way for the next generation of Cardiac Science AEDs.

To learn more about the best choice in AEDs, visit our website [www.cardiacscience.com](http://www.cardiacscience.com), or email [aed@cardiacscience.com](mailto:aed@cardiacscience.com), or call **1.866.846.6628** (USA and Canada only).



CARDIAC SCIENCE

## DEFIBRILLATOR

Operation	Semi-automatic
Waveform	Biphasic truncated exponential
Energy (J) range available	Variable Energy (VE) 140J-250J*
Number of protocols	1
Factory default (Nominal)	200VE, 200VE, 200VE*
Shock time cycle	Less than 55 seconds for three shock series
Controls	One-button operation
Voice prompts	Comprehensive voice instructions guide user through rescue process
Visible indicators	SmartGauge <sup>®</sup> Battery Status Indicator, Service Indicator, Electrodes Indicator
Audible alerts	Voice Prompt, System Alert
Synchronized shock	Built in automatic synchronization
Pacemaker pulse detection	Yes
Programmable	No
Pediatric capability	Yes

## ELECTRODES

Minimum combined surface area	228cm <sup>2</sup>
Extended length of lead wire	1.3m
Supplied	Pre-connected to the AED, test for presence
Type	Adult, pre-gelled, self-adhesive, disposable, non-polarized (identical pads can be placed in either position) defibrillation electrodes
Shelf life	2 years

\*USA only.

## BATTERY

Type	IntelliSense <sup>®</sup> lithium
Warranty	1 year, full replacement

## AUTOMATED SELF-TESTS

Daily	Battery, electrode presence, internal electronics, no energy charge, SHOCK/CONTINUE button, and software
Weekly	Battery, electrode presence, internal electronics, partial energy charge, SHOCK/CONTINUE button, and software
Monthly	Battery, electrode presence, internal electronics, full energy charge cycle, SHOCK/CONTINUE button, and software

## EVENT DOCUMENTATION

Type	Internal memory
Internal memory	34 minutes ECG data with event annotation
Playback	Viewable via RescueLink <sup>®</sup> software via PC
Communications	Serial port or USB (via adapter) for PC with Windows
Clock synchronization	Precise synchronization of AED time to dispatch time

## PHYSICAL DIMENSIONS

Height	8 cm (3.3 in)
Width	27 cm (10.6 in)
Depth	31 cm (12.4 in)
Weight	3.10 kg (6.6 lb)

# Ordering Information

## MODEL

9300C-201

## DESCRIPTION

FirstSave AED G3 Automated External Defibrillator

Each AED package includes: (1) defibrillator, (1) IntelliSense lithium battery (9143), (1) pair of defibrillation electrodes, and (1) Quick Start Tool Kit: includes Quick Start Guide, CD-ROM with AED Manual, Training Video, RescueLink, and serial communication cable.

**TO ORDER CALL:**

"World's Largest AED Source"

**AEDSuperstore<sup>®</sup>**  
an Allied 100<sup>™</sup> company

**800.544.0048**  
**help@aeds.com**

**www.AEDSuperstore.com**

1800 US Hwy 51 N • Woodruff, WI 54568 • 715.358.2329



CARDIAC SCIENCE

