

# Carry Case Set-Up Guide for Philips FRx AED

Front View

- 1. Clear pocket for "Call 911" card



Top Compartment

- 2. Clear pocket for Reference Guide
- 3. Compartment for FRx defibrillator



Middle Compartment  
*Lift white tab to open*

- 4. Mesh pouch for infant/child key
- 5. Compartment for preconnected electrode pads



Bottom Compartment  
*Lift black tab to open*

- 6. Compartment for spare battery
- 7. Compartment for spare electrode pads



Back View

- 8. Clear pocket for emergency contact information and/or owner identification
- 9. Mesh port: push here to release AED from compartment

