

## LIFE® O<sub>2</sub> SoftPac™

### EMERGENCY OXYGEN

can be the difference between  
LIFE OR DEATH

- **Light Weight.....Only 6 Pounds!**  
4 1/2" D x 7 1/2" W x 15" H
- **Comfortable Handle & Shoulder Strap**
- **Wall-Mount Hook for easy access**
- **Easy reading constant supply gauge always viewable through window**  
(Other's units must be turned on to activate supply gauge)
- **40 Minute Supply @ 6 Liters Per Minute**  
240 Liter Aluminum "C" Cylinder
- **Refillable.....Low Cost**  
Universal CGA 870 Pin-Index Valve
- **Simple external 2-step instructions**
  - **Inhalator:** administer supplemental oxygen for a breathing victim.
  - **Resuscitator:** oxygen enriched mouth-to-mask CPR for a non-breathing victim.
- **LIFE® CPR Mask** invertably fits both ADULT and CHILD (patented)  
One-way valve and filter provide hygienic protection which eliminates transmission of vomitus and contagious diseases.
- **Crimp-Proof O<sub>2</sub> Hose**
- **SoftPac™ color choice:** white or green
- **5 Year Warranty**

MODEL	FEATURES
LIFE-2-6FF	6 LPM ( <i>Liters Per Minute</i> ) Fixed Flow
LIFE-2-612	6 & 12 LPM - "NORM" & "HIGH" (2 settings) provides up to recommended 100% inspired O <sub>2</sub>
LIFE-2-025	0 to 25 LPM Variable Flow (for EMT's) 12 flow settings - 0, .25, .5, 1, 2, 3, 4, 6, 8, 10, 15, 25

LIFE® O<sub>2</sub> SoftPac™ meets or exceeds all government and FDA requirements:  
a non-prescription device intended for emergency use.



#### TO ORDER CALL:

"World's Largest AED Source"

**AEDSuperstore®**  
an Allied 100™ company

**800.544.0048**  
help@aeds.com



**www.AEDSuperstore.com**

1800 US Hwy 51 N • Woodruff, WI 54568 • 715.358.2329