

ZOLL

AED+PLUS

AED Plus Accessories and Training Products

AED Plus Trainer

Allows students to practice realistic device operation using AHA recommended training scenarios. An identical graphical operator interface and complete audio and visual prompts mimic the AED Plus in its look and operation. Screen messages, audible tones and voice prompts guide the user through all aspects of the rescue from arrival on scene and calling for help, to administering CPR and shock delivery.

The convenient handheld remote control for instructor use is preprogrammed with the eight American Heart Association Heartsaver AED training scenarios and offers a manual scenario function, plus "Attach Pads," "Low Battery," and "Call for Service" simulations. The Training Unit is powered by six C cell batteries or an AC adapter.

CPR-D Training padz

ZOLL's unique CPR-D Training padz and cables allow users to simulate electrode application and hand placement for CPR. When used for CPR manikin practice, the AED Plus Trainer and CPR-D Training padz help train users to perform CPR with the correct compression depth and rate.



AED Plus Trainer and CPR-D Training padz

TO ORDER CALL:

"World's Largest AED Source"

AEDSuperstore[®]
an Allied 100[™] company

800.544.0048
help@aeds.com



www.AEDSuperstore.com

1800 US Hwy 51 N • Woodruff, WI 54568 • 715.358.2329

ZOLL
It's about time.