

ZOLL AED PLUS MODEL COMPARISON

To Order Call:
AEDSuperstore®
World's Largest Automated External Defibrillator Source
800.544.0048
help@aeds.com
www.AEDSuperstore.com
1800 US Hwy 51 N | Woodruff, WI 54568



ZOLL AED PLUS PS MODEL

Available from AED Superstore

Seven (7) Year Factory Warranty

Features a clear LCD Screen showing voice prompts messages, device advisory messages, elapsed time, shock count, and chest compression graph.

Carry case with rear compartment for spare accessories - Protect your investment!



Operating Temperature:
32 to 122 F (0 to 50 C)

Storage Temperature:
- 22 to 158 F (- 30 to 70 C)

Altitude:
-300 to 15,000 ft. (-91m to 4573m)

ECG/EKG Waveform Display Option

Can be used in loud environments or by a rescuer with hearing impairments.



ZOLL AED PLUS PA MODEL

One (1) Year Factory Warranty

No LCD Screen

No Carry Case

Operating Temperature:
50 to 104 F (10 to 40 C)

Storage Temperature:
32 to 122 F (0 to 50 C)

Altitude:
-300 to 7,500 ft. (-91m to 2287m)

No ECG/EKG Display Option

Not optimal for use in loud environments or by a rescuer with hearing impairments.